



## Canespro Fungal Nail Treatment Set

### Uses

Removes the infected parts of the nail in 2-3 weeks.

### Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.  
Step 2: Cover the treated nail with the plaster provided and keep covered for 24 hours.  
Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.  
Repeat steps 1, 2 and 3 daily for 2-3 weeks.

### Format and Ingredients

Treatment set containing:  
10g tube Canespro Urea 40% ointment, 22 waterproof plasters, 1 plastic nail scraper,

Category:

Pack size: 10g

69

us: MD



Communications International Group Ltd 2025,  
Linen Hall, 162-168 Regent Street, London W1B 5TB