



Carnation Corn Caps

Uses

Removal of hard corns. No more than 3 corns to be treated at any one time.

Dosage

Adults over 15 years: a medicated dressing is placed over the corn and left in place. Change the corn cap every 2 days until the corn is easily removed. A maximum of 5 corn caps should be used per corn. Not to be used by: diabetics, pregnant women, anyone allergic to peanut or soya, people with severe circulatory disorders or those under 15 years of age except on medical advice. Contains peanut oil.

Format and Ingredients

Adhesive dressing with felt ring containing paste of:
Salicylic Acid 40% w/w,

Category:
Footcare

Pack size: 5, 10

RRP: £2.19, £3.89

us: (GSL), (GSL)

A division of
**GROUPE
EUROCOM**



Communications International Group Ltd 2025,
Linen Hall, 162-168 Regent Street, London W1B 5TB