


Canespro Fungal Nail Treatment Set



Category
Footcare
Manufacturer
BAYER PLC
Pack size
10g
RRP
£38.29
Legal Status
MD



Uses

Removes the infected parts of the nail in 2-3 weeks.



Format and Ingredients

Treatment set containing:
10g tube Canespro Urea 40% ointment, 22 waterproof plasters, 1 plastic nail scraper,



Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.
Step 2: Cover the treated nail with one of the plasters provided and keep covered for 24 hours.
Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.
Repeat steps 1, 2 and 3 daily for 2-3 weeks.