



Canespro Fungal Nail Treatment Set

Uses

Removes the infected parts of the nail in 2-3 weeks. As well as Canespro Urea 40% ointment, the set includes 22 single use waterproof plasters and 1 plastic nail scraper.

Dosage

Step 1: Apply a thin layer of the ointment to the infected nail.
Step 2: Cover the treated nail with one of the plasters provided and keep covered for 24 hours.
Step 3: Remove the plaster and soak the nail in warm water for approximately 10 minutes. Gently remove the softened parts of the nail using the plastic scraper.
Repeat steps 1, 2 and 3 daily for 2-3 weeks.

Format and Ingredients

Ointment containing:
Urea 40%, Cera Alba, Lanolin, Petrolatum, Silica,

Category:
Footcare, Skin Care and Antiseptics

Manufacturer:
BAYER PLC

Pack size: Set

RRP: £38.29

Legal Status: MD