



## Sterimar Congestion Relief

### Uses

Helps to rapidly decongest the nose in the event of colds and sinusitis by washing out nasal cavities. Washes away thick mucus.

### Dosage

Adults and children from 3 years: prime the spray by briefly pressing the nozzle. From the onset of colds and sinusitis: Keep the head straight and put the nozzle in the nostril, spraying once or twice into each nostril. Let any excess solution flow out and then blow the nose. For a deeper cleansing of the nasal cavities and blocked noses: Keep the head tilted to one side and put the nozzle in the upper nostril, spraying once or twice for 2-3 seconds. Let any excess solution flow out and then blow the nose. Use 2-6 times a day until the nose is no longer blocked or as long as recommended.

### Format and Ingredients

Nasal spray containing:  
100% natural sea water based solution (hypertonic solution enriched with copper),



Communications International Group Ltd 2025,  
Linen Hall, 162-168 Regent Street, London W1B 5TB

50ml, 100ml

9, £7.99

us: MD, MD